



MATERIAL TRANSFER SYSTEMS

FEATURES AND OPTIONS
CONFIGURATOR

CUSTOMER INFORMATION:

COMPANY: _____
CONTACT: _____
ADDRESS: _____

PHONE: _____
CELL: _____
E-MAIL: _____

ISI REPRESENTATIVE & CONTACT INFO:

NAME: _____
COMPANY: _____
PHONE: _____
CELL: _____
E-MAIL: _____

END LOCATION INFORMATION:

◇ Same as Customer Information

COMPANY: _____
CONTACT: _____
ADDRESS: _____

PHONE: _____
CELL: _____
E-MAIL: _____

PROJECT TIMEFRAME:

- ◇ IMMEDIATE NEED (WITHIN 6 MONTHS)
- ◇ CURRENT NEED (6 MONTHS TO ONE YEAR)
- ◇ FUTURE NEED (MORE THAN ONE YEAR)
- ◇ BUDGETARY
- ◇ OTHER, PLEASE SPECIFY: _____

PROJECT INFORMATION:

STANDARD FEATURES:

- ◆ 304L STAINLESS STEEL CONSTRUCTION WITH #4 FINISH
- ◆ 304L STAINLESS STEEL FRAMED DOORS
- ◆ COVED CORNER INTERIOR TRANSITIONS
- ◆ POLYCARBONATE GLASS
- ◆ MECHANICAL INTERLOCK SYSTEM
- ◆ HEAVY DUTY DOOR HARDWARE
- ◆ INTEGRAL BLOWER A/R
- ◆ HEPA FILTRATION A/R
- ◆ MOUNTING HARDWARE A/R
- ◆ FDA APPROVED SILICONE SEALING GASKET



MATERIAL TRANSFER SYSTEMS FEATURES AND OPTIONS CONFIGURATOR

URS AVAILABLE?:

- ◇ YES

MATERIAL OF CONSTRUCTION:

- ◇ 304 STAINLESS STEEL #4 FINISH
- ◇ 316L STAINLESS STEEL #4 FINISH
- ◇ SPECIAL FINISH SPECIFY _____

DESIGN CONFIGURATIONS:

- ◇ PASS THRU AIR LOCKS WITH MAGNETIC INTERLOCKS
- ◇ LAMINAR DOWN FLOW PASS THRU AIR LOCKS
- ◇ CROSS FLOW PASS THRU AIR LOCKS
- ◇ BLOW OFF PASS THRU AIR LOCKS
- ◇ CART PASS THRU
- ◇ CART AIR SHOWERS
- ◇ DECONTAMINATION CHAMBERS
- ◇ PORTABLE ISO-5 TRANSFER CARTS
- ◇ ALPHA / BETA RAPID TRANSFER PORTS
- ◇ CUSTOMER DESIGNED SYSTEMS

DOORS:

- ◇ STAINLESS STEEL FRAMED DOORS
- ◇ SAFETY GLASS DOORS
- ◇ COUNTER-BALANCED VERTICAL SLIDING DOORS
- ◇ AUTO DOORS

SPECIAL FEATURES:

- ◇ CART TRACK AND WALL BUMPERS
- ◇ SLOPED TOPS
- ◇ SPECIALTY FINISHES

INTERLOCK SYSTEMS:

- ◇ ELECTRO-MECHANICAL INTERLOCK
- ◇ ELECTRO-MAGNETIC INTERLOCK
- ◇ KEYED LOCK OUT

EQUIPMENT INSTALLATION**FILTRATION:**

- ◇ RE-CIRCULATING DOWNFLOW HEPA FILTRATION
- ◇ SINGLE PASS DOWNFLOW HEPA FILTRATION

FILTRATION:

- ◇ HEPA 99.99% eff @ .3 microns
 - ◇ ePTFE HYDROPHOBIC HEPA
- ◇ ULPA 99.9995% eff @ .12 microns
 - ◇ ePTFE HYDROPHOBIC ULPA
- ◇ ASHRAE
 - ◇ SPECIFY RATING _____
- ◇ CARBON FILTER
 - ◇ SPECIFY RATING _____

INTEGRAL BLOWER SYSTEM:

- ◇ SUPPLY
- ◇ EXHAUST

INTERNAL LIGHTING:

- ◇ FLUORESCENT
- ◇ INCANDESCENT
- ◇ ULTRA VIOLET
- ◇ LED
- ◇ VAPOR-TIGHT LED
- ◇ EXPLOSION PROOF CLASS 1 DIV 1 OR DIV 2

DRAINAGE:

- ◇ SPILL DRAIN

DRUM INTRODUCTION LIFT:

- ◇ MANUAL HAND CRANK
- ◇ AUTOMATIC ADJUSTMENT
 - ◇ FOOT PEDAL ACTIVATED
 - ◇ PUSH BUTTON ACTIVATED

FACTORY ACCEPTANCE TESTING (FAT)

- ◇ SITE ACCEPTANCE TESTING (SAT)

CLARIFICATION NOTES:



Isolation Systems Inc.