

CONTAINMENT BOOTH FEATURES AND OPTIONS CONFIGURATOR

CUSTOM SIZE OR CONFIGURATION, PLEASE SPECIFY:

MATERIAL OF CONSTRUCTION:

- ◇ 304 STAINLESS STEEL #4 FINISH
- ◇ 316L STAINLESS STEEL #4 FINISH
- ◇ SPECIAL FINISH SPECIFY _____

WALL PANELS:

- ◇ FULL LIGHT CLEAR GLASS WALL PANELS
- ◇ ANTISTATIC CLEAR ACRYLIC WALL PANELS
- ◇ REMOVEABLE CLEAR FLEXIBLE PVC WALLS
- ◇ MODULAR SOFTWALL DESIGNS W/ MODULAR EXTRACTION MODULE

NEGATIVE PRESSURE OPERATION

- ◇ RE-CIRCULATING
- ◇ ONCE THRU / SINGLE PASS
- ◇ FUME EXTRACTION
- ◇ SECONDARY BARRIER SYSTEMS
- ◇ PORTABLE SYSTEMS
- ◇ SMALL BATCH

STANDARD FEATURES:

- ◇ 304L STAINLESS STEEL
- ◇ PRE-WIRED MODULAR ASSEMBLY
- ◇ 4" THICK DOUBLE PAN HOLLOW WALL SECTIONS
- ◇ SERVICE KNOCKOUTS IN WALLS A/R
- ◇ DAMPERED DISCHARGE GRILLE
- ◇ FLUSH MOUNTED LED LIGHTING
- ◇ NON-SPARKING BACKWARD CURVED ALUMINUM IMPELLARS
- ◇ VARIABLE SPEED CONTROL
- ◇ CG SCREEN DOWN FLOW AIR DIFFUSERS
- ◇ MERV-8 EXTRACTION PRE-FILTERS
- ◇ HIGH VOLUME HEPA EXTRACTION FILTERS
- ◇ FILTERS, BLOWERS, LIGHTS & CONTROLS ARE ROOM SIDE ACCESSIBLE
- ◇ VISUAL OPERATIONAL PRESSURE GAUGES
- ◇ LIGHTED WET LOC. RATED ON / OFF SWITCHES
- ◇ MUSHROOM HEAD EMERGENCY STOP SWITCH
- ◇ SINGLE POINT ELECTRICAL CONNECTION

ARTICULATING SECONDARY SHIELD:

- ◇ ACRYLIC
- ◇ POLYCARBONATE
- ◇ XYZ AXIS

SECONDARY SHIELD GLOVE PORT SIZE:

- ◇ 8" CIRCULAR
- ◇ 10" CIRCULAR
- ◇ 12" CIRCULAR

ARTICULATING SECONDARY SHIELD GLOVES:

- ◇ NATURAL RUBBER
- ◇ NEOPRENE
- ◇ HYPALON (CSM)
- ◇ BUTYL
- ◇ LATEX
- ◇ VITON POLYURETHANE

DOORS:

- ◇ SOLID SS DOORS
- ◇ AUTO ROLL UP DOORS
- ◇ STRIP DOORS

SPECIAL ACCESS DOORS PLEASE SPECIFY:

DRUM INTRODUCTION LIFT:

- ◇ MANUAL HAND CRANK
- ◇ AUTOMATIC ADJUSTMENT, FOOT PEDAL OR PUSH BUTTON ACTIVATED

FILTRATION:

- ◇ HEPA 99.99% eff @ .3 microns
- ◇ ULPA 99.9995% eff @ .12 microns
- ◇ ASHRAE eff. PLEASE SPECIFY
- ◇ CARBON FILTER PLEASE SPECIFY RATING
- ◇ ePTFE HYDROPHOBIC HEPA FILTRATION
- ◇ BAG IN / BAG OUT (BIBO)

BLOWER SYSTEM OPTIONS:

- ◇ EXPLOSION PROOF CLASS 1 DIV 1
- ◇ PRESSURE ACTIVATED BLOWER SPEED MODULATION
- ◇ VARIABLE FREQUENCY DRIVE



Isolation Systems Inc.

CONFIGURATOR CONTINUED

INTEGRAL LIGHTING:

- ◇ HIGH INTENSITY BIAx
- ◇ FLUORESCENT
- ◇ INCANDESCENT
- ◇ ULTRA VIOLET
- ◇ VAPOR-TIGHT LED
- ◇ EXPLOSION PROOF CLASS 1 DIV 1 or DIV 2

ELECTRICAL SERVICES:

- ◇ 120 / 220V 60HZ / 1PH
- ◇ 230 / 460V 60HZ / 3PH

ELECTRICAL OUTLETS:

- ◇ GFI DUPLEX
- ◇ EXPLOSION PROOF
- ◇ SEALED / STAINLESS STEEL COVER
- ◇ STRIP OUTLET
- ◇ SPECIAL VOLTAGE _____
- ◇ SPECIFY AMPERAGE _____

CONTROLS:

- ◇ MAGNEHELIC GAGE
- ◇ MINIHELIC GAGE
- ◇ PHOTOHELIC GAGE
- ◇ DIGIHELIC
- ◇ PLC & ALARM SYSTEM
- ◇ REMOTE LOCATED CONTROL PANELS
- ◇ AUDIBLE VISUAL ALARM SYSTEM
- ◇ BMS OPERATION

◇ **DAMPERED DISCHARGE DUCT CONNECTION**

◇ **EXTRACTION PLENUM CIP SYSTEM**

◇ **CHILLED WATER COIL**

SERVICE PENETRATIONS AND FITTINGS:

- ◇ TRI-CLAMP AND CAP _____
- ◇ BULK HEAD FITTINGS _____
- ◇ CABLE PENETRATIONS _____
- ◇ VALVES & PET COCKS _____
- ◇ COMPUTER CONNECTIONS _____
- ◇ QUICK CONNECTS _____
- ◇ FAUCETS _____
- ◇ SPRINKLER PENETRATIONS _____

◇ **INTEGRAL WORK SURFACES**

◇ **SHELVING** _____

◇ **HOOKS** _____

◇ **FACTORY ACCEPTANCE TESTING (FAT)**

◇ **SITE ACCEPTANCE TESTING (SAT)**

◇ **EQUIPMENT INSTALLATION**

CLARIFICATION NOTES OR DRAWING



Isolation Systems Inc.

