

IMPLANT DRUG ASSEMBLY ISOLATOR



ENGINEERED FOR POSITIVE OR NEGATIVE PRESSURE APPLICATIONS

AVAILABLE FEATURES:

- Fully welded 10 ga. 304 or 316L stainless steel construction with #4 finish on both exterior and interior surfaces. All welds are ground smooth and polished to matching finish.
- All interior corner transitions are coved for ease of cleaning.
- Pass thru antechambers.
- Support frame is 14 ga 2" sq. 304 stainless steel tube with a brushed #4 finish with leg levelers or casters.
- 3/8" thick clear polycarbonate glass view panel.
- Angled ergonomic view panel.
- Hinged service and equipment introduction lift up front sash with gas cylinder assist.
- 12" dia. round or oval Stainless Steel glove ports with compression seal ring.
- One piece HYPALON rubber gloves and sleeves. Optional materials available for gloves.
- 10" Sunken tub for ergonomic internal equipment placement.
- HEPA filtered 99.99 % eff @ .3 microns. Gel sealed to chamber.
- Butterfly valves for chamber seal during filter change out.
- Lighting is LED type.
- Low pressure audible visual alarm.
- ULPA filtered 99.9995 eff @ .12 microns. Gel sealed to chamber.
- Hydrophobic HEPA & ULPA Filtration available.
- Butterfly valves for chamber seal during filter change out.
- Lighting is LED type.
- Low pressure audible visual alarm.
- Magnehelic gage to monitor HEPA filter pressure differential.
- Magnehelic gage to monitor chamber pressurization.
- Optional photohelic pressure switches / gages.
- Wash down blower and light switches.
- Tri-clamp and cap penetrations into main chamber for service connections A/R.
- Hydraulic auto height adjustment.
- Magnetic Interlocks and safety switches.
- Gas Monitors and Alarm systems.
- Option on board Glove Integrity Tester
- Humidity and Temperature Monitors.
- CI-1 Div-1, CI-1 Div-2, CI-2 Div 1 available.

MODEL	DESCRIPTION	ELECTRICAL	INTERIOR DIMENSIONS
ISO-NPIA-2	Single Station Negative Pressure	120V / 60HZ / 1 PH	42"w x 24"d x 38"h
ISO-NPIA-4	Double Station Negative Pressure	120V / 60HZ / 1 PH	82"w x 24"d x 35"h

**All information herein is subject to change without notice.*

**All catalog dimensions are nominal.*



© 7/2021