

BENCH TOP HEPA FILTERED VENTILATED EXHAUSTING ENCLOSURE



Custom configuration shown with optional equipment assistance jib lift and BIBO HEPA Filter housing.

STANDARD FEATURES:

- 304 stainless steel construction with #4 finish.
- Clear Anti-static acrylic enclosure panels.
- Backward curved low profile motorized impellers with non sparking aluminum impeller.
- HEPA filtration 99.99% efficient @ .3 microns.
- Pre-filtration Merv-8 eff Ashrae pleated disposable.
- Integral flush mounted LED lighting.
- Face Velocity airflow monitor with integral alarm.
- Minihelic Differential Pressure Gauge to monitor filter Static Pressure.
- Wet location rated on / off switches for power and lighting
- 12' power cord with hospital rated plug.
- 120 V 60 hz 1 ph operation
- Integral lifting Lugs

NOTE:

- All filters are PSL challenged at factory to insure media integrity.
- Unit performance meets or exceeds standards as specified in Fed. Std. 209E / ISO-5 and IES-RP-CC-002.4.
- All electrical wiring is in conformance with NEMA and NEC standards.

OPTIONAL FEATURES:

- 316L stainless steel construction with #4 finish.
- Stainless steel frame.
- Hydrophobic HEPA filtration.
- ULPA filtration 99.9995% efficient @ .12 microns.
- Optional voltages in 230v-60hz-1ph.
- Iris port entry penetrations.
- Magnehelic Differential Pressure Gauge.
- Audio / visual alarms.
- Bag In Bag Out (BIBO) HEPA filter housing.
- Optional equipment assistance jib lift—push button height adjustment and 180° jib swing.
- Custom sizes and configurations.

MODEL	WORK AREA	OA DIMENSIONS
VEE-315	36"d x 60"w x 42"h	38"d x 61"w x 73"h
VEE-322	36"d x 90"w x 42"h	38"d x 91"w x 73"h
VEE-330	36"d x 120"w x 42"h	38"d x 121"w x 73"h

**All information herein is subject to change without notice.*

**All catalog dimensions are nominal.*



© 8/2022