

ISOLATOR GLOVE TESTER



STANDARD FEATURES:

- 304 Stainless Steel construction.
- Push button Start / Stop.
- Photohelic Differential pressure Switch.
- Adjustable settings.
- Integral compressor. No compressed air required.
- Illuminated Pass / Fail indicator lights.

OPTIONAL FEATURES:

- 316L Stainless Steel construction.
- Integrated HMI controls.

MODEL - ROUND PORT	MODEL—OVAL PORT	SIZE
ISO-GTR-8	ISO-GTO-8	8" GLOVE PORT
ISO-GTR-10	ISO-GTO-10	10" GLOVE PORT
ISO-GTR-12	ISO-GTO-12	12" GLOVE PORT

**All information herein is subject to change without notice.*



© 7/2021